

Safety Data Sheet SUPERB Universal Adhesive

Section 1. Identification of Substance/Preparation and Company/Undertaking

Product Name: SUPERB Universal Adhesive
Manufacturer: Apex Dental Materials, Incorporated
330 Telser Road
Lake Zurich, IL 60047
Emergency Telephone Number: 1-877-273-9123 or 1-847-719-1133
Telephone Number for Information: 1-877-273-9123 or 1-847-719-1133
Date Prepared: October 2nd, 2013

Section 2. Hazards Identification

Irritates the eyes. May cause allergic reaction on skin, flammable

Section 3. Composition/Information on Ingredients

Ingredient	CAS No.	OSHA PEL	ACGIH TWA	Other
Ethanol	64-17-5	1000ppm	N/D	N/A
Acetone	67-64-1	1000ppm	N/D	N/A
Hydroxyethyl methacrylate	868-77-9	N/D	N/D	
2-Methyl-2-Propenoic Acid (1 Methyleneethylidene) Bis (4,1-Phenyleneoxy (2-Hydroxy-3, 1-propanediyl)) Ester	1565-94-2	N/A	N/A	N/A

Section 4. First-aid Measures

Inhalation: Provide fresh air and rest
Skin Contact: Flush with copious amounts of water
Eye Contact: Flush with running water for a minimum of 15 minutes. Seek medical attention
Ingestion: Drink 2 or more glasses of water or milk, seek medical attention. Do not induce vomiting

Section 5. Fire-fighting Measures

Fire Fighting Media: Dry chemical powder, foam or carbon dioxide
Fire and Explosion Risks: Upon combustion, irritating gases may form

Section 6. Accidental Release Measures

Precautions: Prevent skin/eye contact
Method for Clean-up: Wipe with damp cloth, discard appropriately

Section 7. Handling and Storage

Use according to directions for use. Store at room temperature in the original container provided by the manufacturer. Avoid unnecessary exposure to the product.

Section 8. Exposure Controls/Personal Protection

Eye Protection: Wear protective eyewear
Hand Protection: Wear protective gloves
Respiratory Protection: Avoid prolonged breathing of vapors
Skin Protection: Avoid contact, if contact does occur rinse with copious amounts of water

Section 9. Physical/Chemical Properties

Appearance and Odor: Pale yellow liquid, acetone / ethanol odor
Boiling Point: > 50° Celsius
Vapor Pressure (mm Hg) N/A
Vapor Density (Air = 1) N/A
Specific Gravity (Water = 1) N/A
Melting Point < 0° Celsius
Evaporation Rate N/A
Solubility in Water Soluble

Section 10. Stability and Reactivity

Stable. Avoid high temperatures

Section 11. Toxicological Information

Eye Contact: Irritates the eyes
Skin Contact: May irritate the skin
Sensitization: May cause allergic reaction upon contact with the skin
Other: N/A

Section 12. Ecological Information

Use as directed. Follow local law governing disposal of chemicals.

Section 13. Disposal Considerations

Follow local law and regulations for disposal of this material. If none exist follow state or regional laws and regulations

Section 14. Transport Information

Follow local law and regulations for disposal of this material. If none exist follow state or regional laws and regulations

Section 15. Regulatory Information

Classification: R11 Highly flammable., R36 Irritating to eyes., R38 Irritating to skin., R43 May cause sensitization by skin contact., F; Flammable., Xi; Irritant
Safety: S26 In case of contact with eyes, rinse immediately with plenty of water and seek medical advice. S45 In the event of an accident call immediately a physician. If possible show him/her this document.

Section 16. Other Information

Revision Date: 10/02/13
Revision Number: 1

NOTE: This information is believed to be correct as of the date of this SDS. The manufacturer makes no representation as to the completeness or accuracy of this information and supplies it on the condition that the person receiving it will make their own determination as to its suitability or their purposes prior to use. In no event will the manufacturer be responsible for damages of any nature whatsoever resulting from the use of or reliance upon the information contained herein.

For technical assistance please call toll-free: 1-877-APEX-123